


Fluke AS200-R Probe Accessory Set, Red




Product overview: Fluke AS200-R Probe Accessory Set, Red

Probe Accessory Set for usage with VPS200 & PM8918 Series Scope Probes.

Specifications: Fluke AS200-R Probe Accessory Set, Red

Specifications	
Max. input voltage	600 V CAT III, 1000V CAT II per IEC 1010-2-031
Max. voltage to ground	600 V CAT III, 1000V CAT II per IEC 1010-2-031
Compatible with	VPS100, VPS101, VPS121, VPS200-Series and PM8918-Series probes

Ordering information



Fluke AS200 R

Set includes 600V CAT III / 1000V CAT II rated:

- Industrial Alligator for Probe
 - 2-mm Test Probe for Probe
 - 2-mm Test Probe for Banana
 - Industrial Alligator for Banana
 - Ground Lead 4-mm Banana
-

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**