

TECHNICAL DATA

# Fluke RS200 Probe-Accessories Replacement Set




## Key features

Compatible with PM8918, VPS100, VPS101, VPS121 and VPS200-series probes.

## Product overview: Fluke RS200 Probe-Accessories Replacement Set

Contains replacement parts for usage with the VPS200 & PM8918 Scope Probes.

## Ordering information



### Fluke RS200

ul>

- Set includes 600V CAT III / 1000V CAT II rated:
  - Hook Clip for Probe Tip (2 Red, 1 Gray)
  - 4-mm Test Probe for Probe (1 Red, 1 Gray)
  - Ground Lead with Hook Clip (1 Red, 1 Gray)
  - Ground Lead with Mini Alligator (2 Black).
  - Ground Spring for Probe Tip (5x) 30V CAT II Rated
- Max. Input Voltage: 600V CAT III, 1000V CAT II per IEC 1010-2-031
- Max. Voltage to Ground: 600V CAT III, 1000V CAT II per IEC 1010-2-031



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**