

TECHNICAL DATA

# Fluke TL970 Hook and Pincer Set



## Key features

Handy all-in-one kit that includes 3 Test Lead Sets:

- TL940 Mini-Hook Test Lead Set
- TL950 Mini-Pincer Test Lead Set
- TL960 Micro-Hook Test Lead Set
- One year warranty

## Ordering information



### Fluke TL970

---

Handy all-in-one kit that includes:

- 3 Test Lead Sets:
    - TL940 Mini-Hook Test Lead Set
    - TL950 Mini-Pincer Test Lead Set
    - TL960 Micro-Hook Test Lead Set
-

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**