

Fluke CPAK Magnetic Hanging Tool Pouch



Key features

- Heavy duty polyester pouch with magnetic hanger
- One main pocket and one secondary front pocket
- 10 slots in main pocket for screwdrivers and pliers

Product overview: Fluke CPAK Magnetic Hanging Tool Pouch

Keep all your tools protected and at reach with the Fluke CPAK Magnetic Hanging Tool Pouch. The CPAK Pouch keeps your hand tools organized and together within your tool bag. No need to search for a pair of pliers that fell to the bottom when they're all bundled together in this pouch. The two-magnet hanging strap attaches the pouch to a steel cabinet door or metal pipe. A zipper closure on the main pocket keeps 10 tool slots organized for screwdrivers, pliers, and small test tools. A secondary front pocket is ideal for safety glasses and ear plugs, with a hook and loop closure to keep them in.

Specifications: Fluke CPAK Magnetic Hanging Tool Pouch

Size	11 in x 13 in 28 cm x 33 cm
One main pocket, one secondary front pocket. 10 slots in main pocket.	
Heavy duty polyester	
Two rare-earth magnets for hanging strap	

Ordering information



CPAK

CPAK magnetic hanging tool pouch

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/en-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**