

## Fluke C33 Soft Carrying Case



### Key features

- Zippered carrying case with durable vinyl exterior
- Inner pocket to store test leads or small accessories
- Belt loop for easy carrying

### Product overview: Fluke C33 Soft Carrying Case

The C33 Soft Carrying Case is a zippered tool case with a durable, high-quality vinyl exterior. It features an inner pocket that's perfect for storing your test leads with the tools or keeping small accessories organized. The exterior of the case has a belt loop so you can keep your tools on you when you're busy. Grab your tool and accessories in one and get to work quickly.

### Recommended for use with:

- Fluke T100 series (T100, T120 and T140) series Voltage and Continuity Testers
- Fluke 323
- Fluke 325
- Fluke 376 FC
- [Fluke T5 Electrical Tester](#)

## Specifications: Fluke C33 Soft Carrying Case

Dimensions:

27.9 cm x 9.2 cm x 5.1 cm

11 in x 3.6 in x 2 in

## Ordering information



### **Fluke C33**

Fluke C33 Soft Carrying Case

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa  
+31 (0)40 267 5100  
In Canada (800)-36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com/en-th](http://www.fluke.com/en-th)

©2021 Fluke Corporation. Specifications subject to  
change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**