

TECHNICAL DATA

Fluke C75 Accessory Case



Key features

- Zippered vinyl carrying case for small accessories
- Ideal for test leads, probes and alligator clips
- Two inner pockets for storage

Product overview: Fluke C75 Accessory Case

Having the right accessories on-hand is important to make measurements quickly and safely. The C75 Accessory Case is a zippered carrying case designed to keep all the small accessories together in one easy-to-find place. No more digging around your tool bag to find a set of alligator clips or that one test lead cap that fell off. The interior of the bag has two pockets inside so the right accessory will always be at hand.

Specifications: Fluke C75 Accessory Case

Dimensions:

- 179 x 103 x 26 mm
- 7 x 4 x 1 in

Ordering information



Fluke C75

Fluke C75 Accessory Case

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/en-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**