

Fluke CAMO-37 Camouflage Carrying Case



Key features

- Zippered carrying case
- Belt loop
- Inner pocket for storage

Product overview: Fluke CAMO-37 Camouflage Carrying Case

Get the right case for your Fluke tools. The CAMO-37 camouflage carrying case is a zippered carrying case with durable high quality polyester exterior. It has a belt loop and an inner pocket to store test leads or small accessories. It is available in a trendy camouflage design with a 1 year warranty.

Specifications: Fluke CAMO-37 Camouflage Carrying Case

Dimensions (H x W x D)

26.5 x 9 x 3 cm

10.5 x 5.3 x 1.2 in

Ordering information



Fluke CAMO-37

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/en-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**