

TECHNICAL DATA

Fluke 700TTH5K Premium Transmitter Test Hose Kit



Product overview: Fluke 700TTH5K Premium Transmitter Test Hose Kit

Easy to use premium pressure transmitter test hose kits for use with a calibrator or test pump

- 5,000 psi rating
- Included 1/8" NPT connection for permanent mount to a calibrator or test pump
- No tools required, 1/4" NPT male connector for connection to a transmitter

Ordering information



Fluke 700TTH5K

Includes:

- One (1) premium hose kit and adapters
- Premium/no tools required 1/4" NPT transmitter connection fitting

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/en-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
10/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**