

TECHNICAL DATA

Fluke L211 Probe Light Kit



Product overview: Fluke L211 Probe Light Kit

Special value combo pairs the new L200 Probe Light with one of the most popular test leads by Fluke and a handy carrying case.

Ordering information



Fluke L211

Fluke L211 Probe Light Kit

L211 Probe Light Kit includes:

- L200 Probe Light, TL71 Premium DMM Test Lead Set and C75 Case
- Small, rugged light easily attaches to any Fluke test probe
- Flexible silicone test leads are heat and cold resistant
- Zippered case has two inside pockets
- Probe light frees both hands for work
- 120 hours of battery life
- Uses 2 3-volt watch batteries

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/en-th

©2022 Fluke Corporation. Specifications subject to
change without notice.
01/2022

**Modification of this document is not permitted
without written permission from Fluke Corporation.**