

# Fluke TLK287 Electronics Master Test Lead Set



## Key features

Precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points. Micro-grabbers and leads to enable fine SMD test. Mini-alligator clips, grabbers, and probes for virtually every electronics need.

## Product overview: Fluke TLK287 Electronics Master Test Lead Set

Fluke TLK287 Electronics Master Test Lead Set is ideal for electronic test and design. The set includes precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points. Micro-grabbers and leads enable fine SMD testing. Mini-alligator clips, grabbers and probes for virtually every electronics need. And modular test leads and lead-probe combinations for versatility and couplers for long lead applications.

## Ordering information



### Fluke-TLK287

---

Includes:

- Precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points
  - Micro-grabbers and leads to enable fine SMD test
  - Mini-alligator clips, grabbers, and probes for virtually every electronics need
-

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa  
+31 (0)40 267 5100  
In Canada (800)-36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com/en-th](http://www.fluke.com/en-th)

©2022 Fluke Corporation. Specifications subject to  
change without notice.  
01/2022

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**