

TECHNICAL DATA

Fluke Networks MT-8203-20 IntelliTone Test Leads



Key features

- Test leads with Bed of Nails for Pro 100/200
- Two Banana Jack test leads with piercing pin and alligator clips
- U-shaped nose for blade-shaped terminals
- Notched jaw and serrated teeth for gripping
- Large piercing needle for large gauge wires

Product overview: Fluke Networks MT-8203-20 IntelliTone Test Leads

IntelliTone™ MT-8203-20 test leads featuring an insulation-piercing clip on one end and a shrouded, 4-mm right-angle banana plug on the other. The clip features a bed-of-nails for piercing 22-28 AWG wires.

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/en-th

©2022 Fluke Corporation. Specifications subject to
change without notice.
01/2022

**Modification of this document is not permitted
without written permission from Fluke Corporation.**