

Fluke Networks Remote ID Kit for MS-POE MicroScanner™ with Identifiers #2-7



Key features

- Remote ID kit for use with the MicroScanner™ PoE Cable Verifier
- Far-end identifiers #2-7
- Test telephone jacks, Ethernet jacks and CATV outlets

Product overview: Fluke Networks Remote ID Kit for MS-POE MicroScanner™ with Identifiers #2-7

Remote ID kit for MicroScanner PoE, includes #2-7 IDs, this is an optional accessory. Locators help identify connections at patch panels and shows the number of the ID connected to the far end of the cabling to determine network layout, test telephone, Ethernet, and CATV outlets.

Ordering information



Remote ID Kit

Fluke Networks Remote ID Kit for MS-POE MicroScanner™ with Identifiers #2-7



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**