

TECHNICAL DATA

Fluke BTL20 Interactive Battery Analyzer Test Probe



Product overview: Fluke BTL20 Interactive Battery Analyzer Test Probe

Intelligent test probe set for Fluke BT520 Battery Analyzer (temperature measurement function not included). The BTL20 has a built in LCD display, tri-light LED for pass/fail warning indication, voice count for battery ID and integrated save button. Comes with both long and short probe extenders for use in hard to reach battery cabinets.

Ordering information



Fluke BTL20

Interactive Battery Analyzer Test Probe

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**