

TECHNICAL DATA

# Fluke BTL21 Interactive Battery Analyzer Test Probe with Temperature Sensor



## Product overview: Fluke BTL21 Interactive Battery Analyzer Test Probe with Temperature Sensor

Intelligent test probe set with integrated IR temperature sensor designed for Fluke BT521 Battery Analyzer. The BTL21 has a built in LCD display, tri-light LED for pass/fail warning indication, voice count for battery ID and integrated save button. Comes with both long and short probe extenders for use in hard to reach battery cabinets.

## Ordering information



### **Fluke BTL21**

Interactive Battery Analyzer Test Probe with Temperature Sensor

---



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
11/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**