


TECHNICAL DATA

Fluke MA190 Accessory Kit for Medical and Video Applications



Key features

This kit includes:

- 50 Ohm BNC insulated feedthrough terminator
- 50 Ohm BNC terminator with 10:1 signal attenuation
- 1 Ohm current shunt for current measurements
- Safety-designed BNC cable, 1.5 m (5 ft)
- Insulated BNC (f) to 4 mm-banana plug adapter
- Dual 4 mm banana receptacles (1 red, 1 black)

Product overview: Fluke MA190 Accessory Kit for Medical and Video Applications

The MA190 accessory set enables interconnection of the Fluke ScopeMeter 190 Series for use in the field of Medical Imaging and Video systems.

Ordering information



Fluke MA190

Accessory kit for medical and video applications

Fluke MA190 Accessory Kit for Medical and Video Applications

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**