

TECHNICAL DATA

# Fluke Networks MS2-MAG-KIT Magnetic Strap Attachment and Spare Holster



## Key features

- Attaches to any magnetic surface using the strong rare-earth magnet
- Built in hole to attach the strap to j-hooks in your data rack
- Folds the strap to create a simple stand for the meter
- Works with MicroScanner™, SimpliFiber® Pro, FI-500, MultiFiber™ Pro

## Product overview: Fluke Networks MS2-MAG-KIT Magnetic Strap Attachment and Spare Holster

Holster with magnetic strap for MicroScanner™ cable tester, FI-500 FiberInspector™ Micro, SimpliFiber® Pro, and MultiFiber™ Pro

## Ordering information



### **MS2-MAG-KIT**

Fluke Networks MS2-MAG-KIT Magnetic Strap Attachment and Spare Holster

---



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**