


TECHNICAL DATA

## Fluke PM9094/101 Mini Testhook Set for PM8918 series probes




### Product overview: Fluke PM9094/101 Mini Testhook Set for PM8918 series probes

Set consists of: 4 Mini Testhooks, 2 red, 2 gray. Max. voltage: 30 V rms.

**Note:** When using the mini Test-hooks any potential must not exceed 42 Vpk/30 Vrms

## Ordering information



### **Fluke PM9094/101**

Mini Testhook Set for PM8918 series probes

---



**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**