

Fluke VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)



Product overview: Fluke VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)

Double insulated voltage probe for ScopeMeter190 series.

Specifications: Fluke VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)

Specifications	
Attenuation	10x
Bandwidth	200 MHz
Max. input voltage	1000 V
Max. voltage to ground	600 V CAT III 1000 V CAT II per IEC / EN61010-031
Input resistance	10 Mohm
Input capacitance	14 pF
Compensation range	10 to 22 pF



Cable length	1.2m
Warranty	1 Year
Specifications apply when the probe is connected to Fluke 190 Series Scopemeter.	

Ordering information



Fluke VPS200 G

Voltage Probe Set, Gray 200 MHz (one gray)

Each probe set includes:

- 600 V CAT III / 1000 V CAT II rated:
 - 10:1 Voltage Probe
 - Hook Clip for Probe Tip, 4-mm Test Probe for Probe
 - Ground Lead with Hook Clip
 - Ground Lead with Mini Alligator
 - The red probe set also includes a 30 V CAT II rated Ground Spring for Probe Tip
-

Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
12/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**