

TECHNICAL DATA

Fluke TL970 Hook and Pincer Set



Key features

Handy all-in-one kit that includes 3 Test Lead Sets:

- TL940 Mini-Hook Test Lead Set
- TL950 Mini-Pincer Test Lead Set
- TL960 Micro-Hook Test Lead Set
- One year warranty

Ordering information



Fluke TL970

Handy all-in-one kit that includes:

- 3 Test Lead Sets:
 - TL940 Mini-Hook Test Lead Set
 - TL950 Mini-Pincer Test Lead Set
 - TL960 Micro-Hook Test Lead Set
-



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**