

Fluke TLK287 Electronics Master Test Lead Set




Key features

Precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points. Micro-grabbers and leads to enable fine SMD test. Mini-alligator clips, grabbers, and probes for virtually every electronics need.

Product overview: Fluke TLK287 Electronics Master Test Lead Set

Fluke TLK287 Electronics Master Test Lead Set is ideal for electronic test and design. The set includes precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points. Micro-grabbers and leads enable fine SMD testing. Mini-alligator clips, grabbers and probes for virtually every electronics need. And modular test leads and lead-probe combinations for versatility and couplers for long lead applications.

Ordering information



Fluke-TLK287

Includes:

- Precision electronic probes with multiple spring-loaded sharp tips to maximize contact with SMD test points
 - Micro-grabbers and leads to enable fine SMD test
 - Mini-alligator clips, grabbers, and probes for virtually every electronics need
-



Fluke. *Keeping your world up and running.®*

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Canada (800) 36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com

©2021 Fluke Corporation.
Specifications subject to change without notice.
11/2021

**Modification of this document is not permitted
without written permission from Fluke Corporation.**