

# VPS220-G - Voltage Probe Set, 100:1, one grey



## Key features

**Each probe set includes:**

- 600 V CAT III / 1000 V CAT II rated 100:1 Voltage Probe
- Hook Clip for Probe Tip
- Ground Lead with Mini Alligator
- Insulation Sleeve
- 30 V CAT II rated Ground Spring for Probe Tip

## Product overview: VPS220-G - Voltage Probe Set, 100:1, one grey

Double insulated voltage probe for ScopeMeter 190 series.


## Specifications: VPS220-G - Voltage Probe Set, 100:1, one grey

Specification	
Attenuation	100:1
Bandwidth	200 MHz
Max. input voltage	1000 V CAT II or 600 V CAT III



Max. voltage to ground	1000 V CAT II or 600 V CAT III per IEC / EN61010-031
Input resistance	100 MΩ
Input capacitance	4.6 pF
Compensation range	10 to 25 pF
Cable length	1.2 m

## Ordering information



### **Fluke VPS220-G**

Voltage Probe Set, 100:1, one grey

---

Each probe set includes:

- 600 V CAT III / 1000 V CAT II rated 100:1 Voltage Probe
  - Hook Clip for Probe Tip
  - Ground Lead with Mini Alligator
  - Insulation Sleeve
  - 30 V CAT II rated Ground Spring for Probe Tip
-

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**