


TECHNICAL DATA

# Fluke Networks CLIP-SET RJ45 To 8-Clip Test Leads



## Key features

- Clip set includes RJ45 to 8-clip test lead
- 8 clips are grouped into 4 sets of 2 clips each by color code
- Compatible with MicroScanner™ and CableIQ

## Product overview: Fluke Networks CLIP-SET RJ45 To 8-Clip Test Leads

Clip set includes RJ-45-to-8-clip test lead for CableIQ™ or MicroScanner

## Ordering information



### **CLIP-SET RJ45**

Fluke Networks CLIP-SET RJ45 To 8-Clip Test Leads

---

Includes:

- RJ-45-to-8-clip test lead for CableIQ™ or MicroScanner
-

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Canada (800) 36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

©2021 Fluke Corporation.  
Specifications subject to change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**