

TECHNICAL DATA

Fluke C33 Soft Carrying Case



Key features

- Zippered carrying case with durable vinyl exterior
- Inner pocket to store test leads or small accessories
- Belt loop for easy carrying

Product overview: Fluke C33 Soft Carrying Case

The C33 Soft Carrying Case is a zippered tool case with a durable, high-quality vinyl exterior. It features an inner pocket that's perfect for storing your test leads with the tools or keeping small accessories organized. The exterior of the case has a belt loop so you can keep your tools on you when you're busy. Grab your tool and accessories in one and get to work quickly.

Recommended for use with:

- Fluke T100 series (T100, T120 and T140) series Voltage and Continuity Testers
- Fluke 323
- Fluke 325
- Fluke 376 FC
- [Fluke T5 Electrical Tester](#)

Specifications: Fluke C33 Soft Carrying Case

Dimensions:

27.9 cm x 9.2 cm x 5.1 cm

11 in x 3.6 in x 2 in

Ordering information



Fluke C33

Fluke C33 Soft Carrying Case

Fluke. *Keeping your world up and running.®*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

Representative office of Fluke South East Asia Pte Ltd

C/O Danaher Vietnam
Green Power Tower, 11th Floor Unit 2
35 Ton Duch Thang Street, District 1
Ho Chi Minh City
Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: info.asean@fluke.com

www.fluke.com/vn-vi

©2021 Fluke Corporation. Specifications subject to change without notice.
12/2021

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

Modification of this document is not permitted without written permission from Fluke Corporation.