

## TECHNICAL DATA

## Fluke C75 Accessory Case



### Key features

- Zippered vinyl carrying case for small accessories
- Ideal for test leads, probes and alligator clips
- Two inner pockets for storage

### Product overview: Fluke C75 Accessory Case

Having the right accessories on-hand is important to make measurements quickly and safely. The C75 Accessory Case is a zippered carrying case designed to keep all the small accessories together in one easy-to-find place. No more digging around your tool bag to find a set of alligator clips or that one test lead cap that fell off. The interior of the bag has two pockets inside so the right accessory will always be at hand.

### Specifications: Fluke C75 Accessory Case

Dimensions:

- 179 x 103 x 26 mm
- 7 x 4 x 1 in

## Ordering information



### **Fluke C75**

Fluke C75 Accessory Case

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**

PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

**Representative office of Fluke South East Asia Pte Ltd**

C/O Danaher Vietnam

Green Power Tower, 11th Floor Unit 2

35 Ton Duch Thang Street, District 1

Ho Chi Minh City

Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: [info.asean@fluke.com](mailto:info.asean@fluke.com)

[www.fluke.com/vn-vi](http://www.fluke.com/vn-vi)

©2021 Fluke Corporation. Specifications subject to change without notice.

12/2021

**Modification of this document is not permitted without written permission from Fluke Corporation.**