

TECHNICAL DATA

# Fluke 700TTH5K Premium Transmitter Test Hose Kit



## Product overview: Fluke 700TTH5K Premium Transmitter Test Hose Kit

Easy to use premium pressure transmitter test hose kits for use with a calibrator or test pump

- 5,000 psi rating
- Included 1/8" NPT connection for permanent mount to a calibrator or test pump
- No tools required, 1/4" NPT male connector for connection to a transmitter

## Ordering information



### **Fluke 700TTH5K**

Includes:

- One (1) premium hose kit and adapters
- Premium/no tools required 1/4" NPT transmitter connection fitting

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**

PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (905) 890-7600

From other countries +1 (425) 446-5500

**Representative office of Fluke South East Asia Pte Ltd**

C/O Danaher Vietnam

Green Power Tower, 11th Floor Unit 2

35 Ton Duch Thang Street, District 1

Ho Chi Minh City

Vietnam

Tel: +84-8-2220-5371 (ext 103)

Email: [info.asean@fluke.com](mailto:info.asean@fluke.com)

[www.fluke.com/vn-vi](http://www.fluke.com/vn-vi)

©2021 Fluke Corporation. Specifications subject to change without notice.

11/2021

**Modification of this document is not permitted without written permission from Fluke Corporation.**