

## TECHNICAL DATA

# Fluke TLS430 – Test Leads and Alligator Clips (4 black, 1 green)



## Key features

- Test lead replacement set for the 430 Series Power Quality Analyzer
- Safety rating: Cat III 1000 V, CAT IV 600 V

## Ordering information



### **Fluke TLS430**

Test Leads and Alligator Clips (4 black, 1 green)

---

Set includes five test leads and five alligator clips (4 black, 1 green).

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa  
+31 (0)40 267 5100  
In Canada (905) 890-7600  
From other countries +1 (425) 446-5500

**Representative office of Fluke South East Asia Pte Ltd**

C/O Danaher Vietnam  
Green Power Tower, 11th Floor Unit 2  
35 Ton Duch Thang Street, District 1  
Ho Chi Minh City  
Vietnam  
Tel: +84-8-2220-5371 (ext 103)  
Email: info.asean@fluke.com  
www.fluke.com/vn-vi  
©2021 Fluke Corporation. Specifications subject to  
change without notice.  
12/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**