


TECHNICAL DATA

Fluke ES-162P3-2 3-Pole Stake Kit



Product overview: Fluke ES-162P3-2 3-Pole Stake Kit

Needed to perform 3-Pole Fall-of-Potential testing (Stake testing). Used with both the Fluke 1623-2 and Fluke 1625-2. This is a subset of what is included in the Fluke 1625-2 Kit.

Ordering information



Fluke ES 162P3-2

- (3) Stakes
- (1) 50 m cable reel of red wire
- (1) 25 m cable reel of blue wire

Fluke. *Keeping your world up and running.®*

Fluke Europe B.V.

P.O. Box 1186
5602 BD Eindhoven
The Netherlands
www.fluke.com/en

©2022 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
01/2022

For more information call:

In Middle East/Africa
+31 (0)40 267 5100

**Modification of this document is not permitted
without written permission from Fluke Corporation.**