

TECHNICAL DATA

Fluke 80PK-22 SureGrip™ Immersion Temperature Probe



Key features

- Type-K thermocouple for use in liquids and gels
- Measurement range: -40 to 1090 °C
- One year warranty

Product overview: Fluke 80PK-22 SureGrip™ Immersion Temperature Probe

Fluke 80PK-22 SureGrip™ Immersion Temperature Probe is designed for use in liquid and gels. It can also be used as a general purpose temperature probe. The thermocouple junction is protected from tip to handle by an Inconel sheath. The 80PK-22 can be used with any temperature-measuring instrument that is designed to accept type K thermocouples and has a miniature connector input.

Ordering information



Fluke 80PK 22

Fluke 80PK-22 SureGrip™ Immersion Temperature Probe

