

# Fluke TLK-225 SureGrip™ Master Accessory Set



## Key features

The Fluke TLK-225 SureGrip™ Master Accessory Set is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes.

Set includes:

- TL224 Test Lead Set
- TP175 TwistGuard™ Test Probes
- AC285 SureGrip™ Alligator Clips
- AC220 SureGrip™ Alligator Clips
- AC280 SureGrip™ Hook Clips
- AC283 SureGrip™ Pincer Clips
- A convenient six-pocket roll up pouch
- CAT III 1000 V, CAT IV 600 V safety rating

## Product overview: Fluke TLK-225 SureGrip™ Master Accessory Set

Test leads are an integral part of the complete measurement system. Worn out test leads can cause inaccurate readings and pose serious shock or electrocution hazard. The Fluke TLK-225 SureGrip™ Master Accessory Set with 13 items is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes in a convenient six pocket, roll up pouch.

## Ordering information



### Fluke TLK225

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Europe B.V.**

P.O. Box 1186  
5602 BD Eindhoven  
The Netherlands  
[www.fluke.com/en](http://www.fluke.com/en)

©2021 Fluke Corporation. All rights reserved.  
Data subject to alteration without notice.  
12/2021

**For more information call:**

In Middle East/Africa  
+31 (0)40 267 5100

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**