

TECHNICAL DATA

Fluke TL910 Electronic Test Probes



Key features

- Very small 1 mm tips to access hard to reach test points
- Probe tip length up to 100 mm, test lead length: 1 m
- 3 sets of spring loaded gold tips and 2 sets of stainless steel tips included
- CAT II 1000 V, 3 A rating
- One year warranty

Ordering information



Fluke TL910

Fluke. *Keeping your world up and running.®*

Fluke Europe B.V.

P.O. Box 1186
5602 BD Eindhoven
The Netherlands
www.fluke.com/en

©2021 Fluke Corporation. All rights reserved.
Data subject to alteration without notice.
12/2021

For more information call:

In Middle East/Africa
+31 (0)40 267 5100

**Modification of this document is not permitted
without written permission from Fluke Corporation.**