

Datos técnicos

Kit de aspiración, modelo CO-205



Características principales

- Stainless steel sampling tube
- Industrial-grade hand operated aspirator to draw flue sample
- Easy to replace particulate filter
- Specially designed nose cap for connection to the Fluke CO-210/220
- 4' (1.2 m) hose length

Descripción general del producto: Kit de aspiración, modelo CO-205

Sample flue gas with your Fluke CO-210 or CO-220 using the CO-205 Aspirator Kit

Draw flue gas samples up to 700°F (370°C) reliably for measurement with the Fluke CO-210 Carbon Monoxide probe or Fluke CO-220 Carbon Monoxide meter.

The CO-205 aspirator accessory kit contains all the components necessary to provide a clean sample for the Fluke family of gas measuring devices.

Modelos



Fluke CO-205

Fluke CO-205 Aspirator Kit

Fluke. *Manteniendo su mundo en marcha.*

Fluke Corporation
Everett, WA 98206 EE.UU.

Para obtener información adicional En EE. UU. (800)
443-5853

En Europa/Medio Oriente/África
+31 (0)40 267 5100

En Canadá (800)-36-FLUKE
www.fluke.com

Latin America
Tel: +1 (425) 446-5500
www.fluke.com/es-gt

©2022 Fluke Corporation. Reservados todos los derechos. Impreso en los Países Bajos. Información sujeta a modificación sin previo aviso.
01/2022

No está permitido modificar este documento sin autorización por escrito de Fluke Corporation.