

Datos técnicos

## Kit de aspiración, modelo CO-205



### Características principales

- Stainless steel sampling tube
- Industrial-grade hand operated aspirator to draw flue sample
- Easy to replace particulate filter
- Specially designed nose cap for connection to the Fluke CO-210/220
- 4' (1.2 m) hose length

### Descripción general del producto: Kit de aspiración, modelo CO-205

#### Sample flue gas with your Fluke CO-210 or CO-220 using the CO-205 Aspirator Kit

Draw flue gas samples up to 700°F (370°C) reliably for measurement with the Fluke CO-210 Carbon Monoxide probe or Fluke CO-220 Carbon Monoxide meter.

The CO-205 aspirator accessory kit contains all the components necessary to provide a clean sample for the Fluke family of gas measuring devices.

## Modelos



### **Fluke CO-205**

Fluke CO-205 Aspirator Kit

---

**Fluke. Manteniendo su mundo en marcha.**

**Fluke Corporation**  
Everett, WA 98206 EE.UU.

**Para obtener información adicional En EE. UU. (800) 443-5853**

**En Europa/Medio Oriente/África**  
**+31 (0)40 267 5100**

**En Canadá (800)-36-FLUKE**  
**www.fluke.com**

**Latin America**  
Tel: +1 (425) 446-5500  
[www.fluke.com/laam](http://www.fluke.com/laam)

©2022 Fluke Corporation. Reservados todos los derechos. Impreso en los Países Bajos. Información sujeta a modificación sin previo aviso.  
01/2022

**No está permitido modificar este documento sin autorización por escrito de Fluke Corporation.**