

Datos técnicos

TLK281 -1, Automotive Test Lead Kit

Descripción general del producto: TLK281 -1, Automotive Test Lead Kit

Complete set of test and measurement accessories designed especially for automotive specialists. Use the sharp stainless steel probes to contact the wire conductor or pinpoint tests on a circuit or relay, and the large jaw alligator clips to get a good ground on the engine or to connect to the battery terminals. An extra long set of alligator clips easily grip to the alternator BAT terminal or the S terminals. The set comes with silicone test leads and a case to keep it organized.

- Tested to 60 V DC

Kit contains:

- **TP81 Insulation Piercing Clip Set**
- **TL224 Suregrip™ Silicone Test Lead Set**
- **TP220-1 Suregrip™ Test Probe Set**
- **AC220 Suregrip™ Alligator Clip Set**
- **AC285 Suregrip™ Large Jaw Alligator Clip Set**
- **Soft Case**

Fluke. *Manteniendo su mundo en marcha.*

Fluke Corporation
Everett, WA 98206 EE.UU.

Para obtener información adicional En EE. UU. (800) 443-5853

En Europa/Medio Oriente/África
+31 (0)40 267 5100

En Canadá (800)-36-FLUKE
www.fluke.com

Latin America
Tel: +1 (425) 446-5500
www.fluke.com/es-us

©2021 Fluke Corporation. Reservados todos los derechos. Impreso en los Países Bajos. Información sujeta a modificación sin previo aviso.
12/2021

No está permitido modificar este documento sin autorización por escrito de Fluke Corporation.