

Datos técnicos

## P5544 Two Gauge Stand



### Descripción general del producto: P5544 Two Gauge Stand

This adapter mounts directly to the test port of the calibrator and allows for the calibration of two instruments at the same time or the connection of a reference test instrument. Maximum working pressure is 10,000 psi (700 bar).

### Especificaciones: P5544 Two Gauge Stand

Specifications	
Body Material	Brass and Stainless steel
O-Rings	Nitrile
Test Port Adapters	1 set, 1/8, 1/4, 3/8 and 1/2 NPT or BSP and adapters supplied with calibrator
Calibrator Test Connection	1/2 BSP male (mating adapter supplied)

## Modelos

**P5544**

P5544 Two Gauge Stand

---

Two Gauge Stand

---

**Fluke.** *Manteniendo su mundo en marcha.*

**Fluke Corporation**  
Everett, WA 98206 EE.UU.

**Para obtener información adicional En EE. UU. (800)**  
**443-5853**

**En Europa/Medio Oriente/África**  
**+31 (0)40 267 5100**

**En Canadá (800)-36-FLUKE**  
**[www.fluke.com](http://www.fluke.com)**

**Latin America**  
Tel: +1 (425) 446-5500  
[www.fluke.com/es-us](http://www.fluke.com/es-us)

©2021 Fluke Corporation. Reservados todos los derechos. Impreso en los Países Bajos. Información sujeta a modificación sin previo aviso.  
11/2021

**No está permitido modificar este documento sin autorización por escrito de Fluke Corporation.**