



## Modèles



### Fluke-TLK289

Industrial Master Test Lead Set

Includes:

- TL224 test lead set, AC285 SureGrip™ alligator clips, AC220 SureGrip™ alligator clips, AC280 SureGrip™ hook clips and TP220 SureGrip™ industrial test probes for dependable contact with a variety of test points
- TPAK ToolPak™ Magnetic hanger for convenient suspension of your DMM from a metal surface
- 80BK-AK Type thermocouple temperature probe for direct temperature measurement with your DMM's temperature measurement functions
- C116 protective soft case for your DMM and all its accessories

*Soyez à la pointe du progrès avec **Fluke**.*

**Fluke (Switzerland) GmbH**  
Industrial Division  
Hardstrasse 20  
CH-8303 Bassersdorf  
Tel: +41 (0) 44 580 7504  
E-mail: roc.switzerland@fluke.com  
www.fluke.com/fr-ch

©2022 Fluke Corporation. Tous droits réservés.  
Informations modifiables sans préavis.  
01/2022

**La modification de ce document est interdite sans  
l'autorisation écrite de Fluke Corporation.**