

Kit d'accessoires complet SureGrip™ Cordons et sondes Fluke TLK-225



Principales fonctions

The Fluke TLK-225 SureGrip™ Master Accessory Set is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes.

Set includes:

- TL224 Test Lead Set
- TP175 TwistGuard™ Test Probes
- AC285 SureGrip™ Alligator Clips
- AC220 SureGrip™ Alligator Clips
- AC280 SureGrip™ Hook Clips
- AC283 SureGrip™ Pincer Clips
- A convenient six-pocket roll up pouch
- CAT III 1000 V, CAT IV 600 V safety rating

Présentation du produit: Kit d'accessoires complet SureGrip™ Cordons et sondes Fluke TLK-225

Les cordons de mesure font partie intégrante du système de mesure complet. Des cordons de mesure usés peuvent donner des mesures imprécises et présenter des risques sérieux de choc ou d'électrocution. Le kit d'accessoires complet SureGrip™ Fluke TLK-225 et ses 13 éléments est le kit idéal de rechange qui offre de la flexibilité et du confort avec tous les cordons et

sondes SureGrip™ dans une sacoche souple à six compartiments.

Modèles



Fluke TLK225

*Soyez à la pointe du progrès avec **Fluke**.*

Fluke Europe B.V.

P.O. Box 1186
5602 BD Eindhoven
The Netherlands
www.fluke.com/fr

©2021 Fluke Corporation. Tous droits réservés.
Informations modifiables sans préavis.
12/2021

En savoir plus:

Middle East/Africa
+31 (0)40 267 5100

**La modification de ce document est interdite sans
l'autorisation écrite de Fluke Corporation.**