

Data teknis

Fluke 700MTH Premium Transmitter Test Hose Kit



Ikhtisar Produk: Fluke 700MTH Premium Transmitter Test Hose Kit

Kit selang pengujian pemancar tekanan premium yang mudah digunakan untuk koneksi ke metrik 1/4" betina/BSP

- Rating 5.000 psi
- Termasuk koneksi NPT 1/8" untuk pemasangan permanen pada kalibrator atau pompa pengujian
- Konektor jantan yang tidak memerlukan alat untuk koneksi ke metrik 1/4" betina/BSP

Model



Fluke 700MTH

Includes:

- One (1) 5,000 psi rated, premium test hose
 - Premium, no tools required male 1/4" metric/BSP connector
-

Fluke. *Memastikan aktivitas Anda terus berjalan dan beroperasi.*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

www.fluke.com

BUT. FLUKE SOUTH EAST ASIA PTE LTD

Menera Satu Sentra Kelapa Gading #06-05

Jl. Bulevar Kelapa Gading Kav. LA# No. 1

Summarecon Kelapa Gading

Jakarta Utara 14240

Indonesia

Tel: +62 21 2938 5922

Email: info.asean@fluke.com

www.fluke.com/id

©2021 Fluke Corporation. Specifications subject to change without notice.

12/2021

Modification of this document is not permitted without written permission from Fluke Corporation.