

Data teknis

# Fluke TLK-225 SureGrip™ Master Accessory Set



## Fitur utama

The Fluke TLK-225 SureGrip™ Master Accessory Set is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes.

Set includes:

- TL224 Test Lead Set
- TP175 TwistGuard™ Test Probes
- AC285 SureGrip™ Alligator Clips
- AC220 SureGrip™ Alligator Clips
- AC280 SureGrip™ Hook Clips
- AC283 SureGrip™ Pincer Clips
- A convenient six-pocket roll up pouch
- CAT III 1000 V, CAT IV 600 V safety rating

## Ikhtisar Produk: Fluke TLK-225 SureGrip™ Master Accessory Set

Test leads are an integral part of the complete measurement system. Worn out test leads can cause inaccurate readings and pose serious shock or electrocution hazard. The Fluke TLK-225 SureGrip™ Master Accessory Set with 13 items is the perfect replacement kit that provides flexibility and comfort with all SureGrip™ leads and probes in a convenient six pocket, roll up pouch.

## Model



**Fluke TLK225**

---

**Fluke.** *Memastikan aktivitas Anda terus berjalan dan beroperasi.*

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 U.S.A.

**For more information call:**  
In the U.S.A. (800) 443-5853  
In Europe/M-East/Africa  
+31 (0)40 267 5100  
In Canada (800)-36-FLUKE  
From other countries +1 (425) 446-5500  
[www.fluke.com](http://www.fluke.com)

**BUT. FLUKE SOUTH EAST ASIA PTE LTD**  
Menera Satu Sentra Kelapa Gading #06-05  
JI. Bulevar Kelapa Gading Kav. LA# No. 1  
Summarecon Kelapa Gading  
Jakarta Utara 14240  
Indonesia  
Tel: +62 21 2938 5922  
Email: [info.asean@fluke.com](mailto:info.asean@fluke.com)  
[www.fluke.com/id](http://www.fluke.com/id)  
©2021 Fluke Corporation. Specifications subject to  
change without notice.  
11/2021

**Modification of this document is not permitted  
without written permission from Fluke Corporation.**