

Data teknis

Fluke T5-1000 Tester Kit



Ikhtisar Produk: Fluke T5-1000 Tester Kit



Alat pengukur tegangan, kontinuitas, dan arus dengan holster

Model



T5-1000 Kit

Fluke T5-1000 Tester Kit

Includes:

- Fluke T5-1000 Voltage, Continuity and Current Tester
- H5 Electrical Tester Holster (for Fluke T3/T5)
- Fluke 1AC II VoltAlert Voltage Detector
- 4 mm detachable probes
- Instruction sheet

Optional accessories

Description

Fluke H5 Electrical Tester Holster

Rugged fabric holster includes flap for lead storage and built-in belt loop.

Fluke CAMO-37 Camouflage Carrying Case

Get the right case for your Fluke tools. The CAMO-37 camouflage carrying case is a zippered carrying case with durable high quality polyester exterior. It has a belt loop and an inner pocket to store test leads or small accessories. It is available in a trendy camouflage design with a 1 year warranty.

Fluke AC220 SureGrip™ Alligator Clips

SureGrip™ accessories are designed to improve steadiness in slippery hands.

Fluke T5 Tester Accessory Starter Kit

Special value starter kit for the electrical professional who already owns a Fluke T5 Electrical Tester or is thinking of buying one.

Fluke. *Memastikan aktivitas Anda terus berjalan dan beroperasi.*

Fluke Corporation

PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:

In the U.S.A. (800) 443-5853

In Europe/M-East/Africa

+31 (0)40 267 5100

In Canada (800)-36-FLUKE

From other countries +1 (425) 446-5500

www.fluke.com

BUT. FLUKE SOUTH EAST ASIA PTE LTD

Menera Satu Sentra Kelapa Gading #06-05

Jl. Bulevar Kelapa Gading Kav. LA# No. 1

Summarecon Kelapa Gading

Jakarta Utara 14240

Indonesia

Tel: +62 21 2938 5922

Email: info.asean@fluke.com

www.fluke.com/id

©2021 Fluke Corporation. Specifications subject to change without notice.

12/2021

Modification of this document is not permitted without written permission from Fluke Corporation.