

Tekniske data

Fluke 700HTPK Hydraulic Test Pressure Kit



Nøkkelfunksjoner

- Easy to use, rugged and reliable construction, capable of generating hydraulic pressure up to 10,000 psi/690 bar
- Combine with one of the Fluke 700G Precision Pressure Gauges for a complete pressure testing and calibration solution
- In addition, combine with the Fluke 700G/TRACK Software to upload over 8,000 pressure measurements logged in the field to a PC

Produktoversikt: Fluke 700HTPK Hydraulic Test Pressure Kit

Generates pressure up to 10,000 psi, 690 BAR. Includes the 700HTP-2 hydraulic hand pump, test hoses and adapters, and protective hard case.

Modeller



Fluke-700HTPK

Sett for hydraulikktesttrykk med presisjonstrykkmåler for testing, hydraulisk håndpumpe, ett slangesett og adaptere – inkluderer hard koffert

Includes:

- 1 hose kit and adapters
 - Hard case
-

Fluke. *Keeping your world up and running.*®

Fluke Norge AS
Postboks 383
1411 Kolbotn
Tlf: 800 18 227
E-mail: cs.no@fluke.com
www.fluke.no

©2021 Fluke Corporation. Med enerett.
Informasjonen kan endres uten varsel. Vi tar
forbehold om trykkfeil.
12/2021

**Endring av dette dokumentet er ikke tillatt uten
skriftlig tillatelse fra Fluke Corporation.**