

Tekniske data

Fluke 700TTH5K Premium Transmitter Test Hose Kit



Produktoversikt: Fluke 700TTH5K Premium Transmitter Test Hose Kit

Brukervennlig, prima testslangepakke for trykktransmittere, for bruk med kalibratorer og testpumper

Modeller



Fluke 700TTH5K

Includes:

- One (1) premium hose kit and adapters
- Premium/no tools required 1/4" NPT transmitter connection fitting

Fluke. *Keeping your world up and running.*®

Fluke Norge AS
Postboks 383
1411 Kolbotn
Tlf: 800 18 227
E-mail: cs.no@fluke.com
www.fluke.no

©2021 Fluke Corporation. Med enerett.
Informasjonen kan endres uten varsel. Vi tar
forbehold om trykkfeil.
12/2021

**Endring av dette dokumentet er ikke tillatt uten
skriftlig tillatelse fra Fluke Corporation.**