


Tekniske data

## Fluke PM9094/101 Mini Testhook Set for PM8918 series probes




### Produktoversikt: Fluke PM9094/101 Mini Testhook Set for PM8918 series probes

Set consists of: 4 Mini Testhooks, 2 red, 2 gray. Max. voltage: 30 V rms.

**Note:** When using the mini Test-hooks any potential must not exceed 42 Vpk/30 Vrms

## Modeller



### **Fluke PM9094/101**

Mini Testhook Set for PM8918 series probes

---

**Fluke.** *Keeping your world up and running.®*

**Fluke Norge AS**  
Postboks 383  
1411 Kolbotn  
Tlf: 800 18 227  
E-mail: [cs.no@fluke.com](mailto:cs.no@fluke.com)  
[www.fluke.no](http://www.fluke.no)

©2021 Fluke Corporation. Med enerett.  
Informasjonen kan endres uten varsel. Vi tar  
forbehold om trykkfeil.  
12/2021

**Endring av dette dokumentet er ikke tillatt uten  
skriftlig tillatelse fra Fluke Corporation.**