

Dane techniczne

Zestaw - tester elektryczny, holster i wskaźnik 1AC II



Najważniejsze cechy

- T5-1000 electrical tester automatically measures volts AC and volts DC with precise digital resolution
- Includes easy, accurate OpenJaw™ current measurement
- Displays resistance to 1000 Ω
- Includes a continuity beeper for quick go/no-go tests
- All inputs are protected to EN61010-1 CAT III 1000 V CAT IV 600 V. UL, CSA
- 1AC II offers Voltbeat™ technology and continuous self-test so you always know it's working
- Upon detection of voltage, the VoltAlert™ tip glows and beeper sounds
- CAT IV 1000 V for added protection

Omówienie produktu: Zestaw - tester elektryczny, holster i wskaźnik 1AC II

Zestaw zawierający tester napięcia, ciągłości i prądu z holsterem. Dodatkowo gratis 1AC.

Modele



Fluke T5-H5-1AC II Kit

This kit includes:

- Fluke T5-1000
 - Fluke 1AC II
 - H5
-

Fluke. *Keeping your world up and running.*®

Fluke Europe B.V.
P.O. Box 1186
5602 BD Eindhoven
The Netherlands
Tel: +31 4 0267 5406
E-mail cee.cs@fluke.com
www.fluke.pl

©2021 Fluke Corporation. Wszelkie prawa
zastrzeżone.
Dane mogą ulec zmianie bez uprzedzenia.
11/2021

**Modyfikacja niniejszego dokumentu bez pisemnej
zgody Fluke Corporation jest zabroniona.**