

Dados técnicos

## Fluke 700PRV-1 Pressure Relief Valve Kit



### Principais recursos

- Adjustable from 725 to 5800 PSI (50 to 200 bar)
- Preset at 2,900 psi/200 bar
- Repeatability to 10% of nominal setting
- Installs directly into the 700HTP-1 pump

### Informações gerais sobre o produto: Fluke 700PRV-1 Pressure Relief Valve Kit

The Fluke 700PRV-1 consists of a relief valve settable from 725 to 5800 PSI (50 to 200 bar) to be used with the 700HTP-1 Hydraulic Test Pump. The 1/8" NPT male thread to fits directly into the Fluke 700HTP-1 pressure pump.

To be used with the 700HTP-1 Hydraulic Test Pump. The 700PRV uses a multi-turn adjustment screw to set preload on internal disc springs. This relief valve will protect the most commonly used Fluke pressure modules and pressure instruments from damage due to accidental over-pressurization.

### Especificações: Fluke 700PRV-1 Pressure Relief Valve Kit



Especificações mecânicas e gerais	
Tamanho	
Peso	
Garantia	
Material	

## Modelos



### **Fluke 700PRV-1**

Fluke 700PRV-1 Pressure Relief Valve Kit

- One adjustable relief valve
- Instruction sheet

**Fluke.** *Keeping your world up and running.*®

**Fluke Corporation**  
PO Box 9090, Everett, WA 98206 EUA

**Para obter mais informações, ligue para os seguintes números:**

Tel: (11) 3530-8901  
Email: [info@fluke.com.br](mailto:info@fluke.com.br)  
[www.fluke.com/pt-br](http://www.fluke.com/pt-br)

**Fluke do Brasil Ltda**  
CENESP Av. Maria Coelho Aguiar,  
215 – Bloco G, 1º andar  
São Paulo/SP – CEP: 05804-900

©2021 Fluke Corporation. Todos os direitos reservados. Os dados fornecidos estão sujeitos a alterações sem aviso prévio.  
12/2021

**É proibido modificar este documento sem permissão escrita da Fluke Corporation.**