

Dados técnicos

## P5544 Two Gauge Stand



### Informações gerais sobre o produto: P5544 Two Gauge Stand

This adapter mounts directly to the test port of the calibrator and allows for the calibration of two instruments at the same time or the connection of a reference test instrument. Maximum working pressure is 10,000 psi (700 bar).

### Especificações: P5544 Two Gauge Stand

Specifications	
Body Material	Brass and Stainless steel
O-Rings	Nitrile
Test Port Adapters	1 set, 1/8, 1/4, 3/8 and 1/2 NPT or BSP and adapters supplied with calibrator
Calibrator Test Connection	1/2 BSP male (mating adapter supplied)

## Modelos



### **P5544**

P5544 Two Gauge Stand

---

Two Gauge Stand

---

**Fluke.** *Keeping your world up and running.*®

**Fluke Corporation**

PO Box 9090, Everett, WA 98206 EUA

**Para obter mais informações, ligue para os seguintes números:**

Tel: (11) 3530-8901

Email: [info@fluke.com.br](mailto:info@fluke.com.br)

[www.fluke.com/pt-br](http://www.fluke.com/pt-br)

**Fluke do Brasil Ltda**

CENESP Av. Maria Coelho Aguiar,

215 – Bloco G, 1º andar

São Paulo/SP – CEP: 05804-900

©2021 Fluke Corporation. Todos os direitos

reservados. Os dados fornecidos estão sujeitos a alterações sem aviso prévio.

11/2021

**É proibido modificar este documento sem permissão escrita da Fluke Corporation.**