

Технические данные

## VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)



### Обзор прибора: VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)

Double insulated voltage probe for ScopeMeter190 series.

### Характеристики: VPS200-G Voltage Probe Set, Gray 200 MHz (one gray)

Specifications	
Attenuation	10x
Bandwidth	200 MHz
Max. input voltage	1000 V
Max. voltage to ground	600 V CAT III 1000 V CAT II per IEC / EN61010-031
Input resistance	10 Mohm
Input capacitance	14 pF



Compensation range	10 to 22 pF
Cable length	1.2m
Warranty	1 Year
Specifications apply when the probe is connected to Fluke 190 Series Scopemeter.	

## Модели



### Fluke VPS200 G

Voltage Probe Set, Gray 200 MHz (one gray)

---

Each probe set includes:

- 600 V CAT III / 1000 V CAT II rated:
    - 10:1 Voltage Probe
    - Hook Clip for Probe Tip, 4-mm Test Probe for Probe
    - Ground Lead with Hook Clip
    - Ground Lead with Mini Alligator
  - The red probe set also includes a 30 V CAT II rated Ground Spring for Probe Tip
-

**Fluke.** *Keeping your world up and running.*®

ООО "Флюк СИИЭС"  
125993, г. Москва, Ленинградский  
проспект д. 37 к. 9 подъезд 4, 1 этаж,  
БЦ «Аэростар»  
Тел: +7 (495) 664-75-12  
Факс: +7 (495) 664-75-12  
e-mail: info@fluke.ru

© Авторское право 2021 Fluke Corporation.  
Авторские права защищены. Данные могут  
быть изменены без уведомления.  
Самые надежные инструменты в мире  
12/2021

**Не разрешается вносить изменения в данный  
документ без письменного согласия компании  
Fluke Corporation.**