

Технические данные

TL970 Набор крючков и пинцетов



Ключевые особенности

Handy all-in-one kit that includes 3 Test Lead Sets:

- TL940 Mini-Hook Test Lead Set
- TL950 Mini-Pincer Test Lead Set
- TL960 Micro-Hook Test Lead Set
- One year warranty

Модели



Fluke TL970

Handy all-in-one kit that includes:

- 3 Test Lead Sets:
 - TL940 Mini-Hook Test Lead Set
 - TL950 Mini-Pincer Test Lead Set
 - TL960 Micro-Hook Test Lead Set
-

Fluke. Keeping your world up and running.®

ООО "Флюк СИИЭС"
125993, г. Москва, Ленинградский
проспект д. 37 к. 9 подъезд 4, 1 этаж,
БЦ «Аэростар»
Тел: +7 (495) 664-75-12
Факс: +7 (495) 664-75-12
e-mail: info@fluke.ru

© Авторское право 2021 Fluke Corporation.
Авторские права защищены. Данные могут
быть изменены без уведомления.
Самые надежные инструменты в мире
12/2021

**Не разрешается вносить изменения в данный
документ без письменного согласия компании
Fluke Corporation.**