

Технические данные

P5544 Two Gauge Stand



Обзор прибора: P5544 Two Gauge Stand

This adapter mounts directly to the test port of the calibrator and allows for the calibration of two instruments at the same time or the connection of a reference test instrument. Maximum working pressure is 10,000 psi (700 bar).

Характеристики: P5544 Two Gauge Stand

| Specifications | |
|----------------------------|---|
| Body Material | Brass and Stainless steel |
| O-Rings | Nitrile |
| Test Port Adapters | 1 set, 1/8, 1/4, 3/8 and 1/2 NPT or BSP and adapters supplied with calibrator |
| Calibrator Test Connection | 1/2 BSP male (mating adapter supplied) |

Модели



P5544

Two Gauge Stand

Two Gauge Stand

Fluke. *Keeping your world up and running.*®

ООО «Флюк СИИЙЭС»
125993, г. Москва, Ленинградский
проспект д. 37 к. 9 подъезд 4, 1 этаж,
БЦ «Аэростар»
Тел: +7 (495) 664-75-12
Факс: +7 (495) 664-75-12
e-mail: info@fluke.ru

© Авторское право 2021 Fluke Corporation.
Авторские права защищены. Данные могут
быть изменены без уведомления.
Самые надежные инструменты в мире
12/2021

**Не разрешается вносить изменения в данный
документ без письменного согласия компании
Fluke Corporation.**