

ข้อมูลด้านเทคนิค

Fluke C33 Soft Carrying Case



คุณลักษณะสำคัญ

- Zippered carrying case with durable vinyl exterior
- Inner pocket to store test leads or small accessories
- Belt loop for easy carrying

ภาพรวมผลิตภัณฑ์: Fluke C33 Soft Carrying Case

- Zippered carrying case with durable vinyl exterior
- Has belt loop and an inner pocket to store test leads or small accessories
- Companion case to:
 - Fluke T100 series (T100, T120 and T140) series Voltage and Continuity Testers
 - Fluke 330 series Clamp Meters
 - Fluke T5 Electrical Tester
- 1 year warranty

ข้อมูลจำเพาะ: Fluke C33 Soft Carrying Case

Dimensions:

27.9 cm x 9.2 cm x 5.1 cm

11 in x 3.6 in x 2 in

รุ่น



Fluke C33

Fluke C33 Soft Carrying Case

Fluke. ให้โลกของคุณ คงอยู่แต่ละก้าวต่อไป

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.