

ข้อมูลด้านเทคนิค

Fluke 700PRV Pressure Relief Valve Kit



คุณลักษณะสำคัญ

Repeatability \pm 10% of nominal setting. Multi-turn adjustment screw to set preload on internal disc springs.

Included:

- The 700PRV contains two relief valves
 - One adjustable from 725-2900 psi/50-200 bar (preset at 1360 psi/94 bar)
 - One adjustable from 2900-5800 psi/200-400 bar (preset at 5450 psi/376 bar)
- Instruction sheet

ภาพรวมผลิตภัณฑ์: Fluke 700PRV Pressure Relief Valve Kit

The Fluke 700PRV consists of two relief valves (1360 and 5450 psi) to be used with the 700HTP Hydraulic Test Pump. These relief valves will protect the most commonly used Fluke pressure modules from damage due to accidental over-pressurization. 1/4 BSP male parallel thread to fit Fluke 700HTP.

ข้อมูลจำเพาะ: Fluke 700PRV Pressure Relief Valve Kit

Mechanical and General Specifications

Size	40 x 22 mm
------	------------



Weight	60 g per relief valve
Warranty	One year
Material	Type 303 stainless steel

รุ่น**Fluke 700PRV**Fluke 700PRV Pressure Relief Valve Kit

Fluke. ให้โลกของคุณ คงอยู่แต่ละก้าวต่อไป

Fluke Corporation
PO Box 9090, Everett, WA 98206 U.S.A.

For more information call:
In the U.S.A. (800) 443-5853
In Europe/M-East/Africa
+31 (0)40 267 5100
In Canada (800)-36-FLUKE
From other countries +1 (425) 446-5500
www.fluke.com/th-th

©2021 Fluke Corporation. Specifications subject to
change without notice.
12/2021

Modification of this document is not permitted without
written permission from Fluke Corporation.